TEMPLOW L

2101280

**ARMARIO CONGELADOR
UPRIGHT FREEZER**

Indice

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

ENGLISH VERSION

* Stacking mark →


Only load the top shelf up to the stacking mark.
This is important so as to ensure that the air can circulate properly and the temperature is even throughout the interior.

1. GENERAL INSTRUCTIONS

- 1) To prevent injury or damage to the unit, the appliance should be unpacked and set up by two people.
- 2) In the event that the appliance is damaged on delivery, contact the supplier immediately before connecting to the mains.
- 3) To guarantee safe operation, ensure that the appliance is set up and connected as described in these operating instructions.
- 4) Disconnect the appliance from the mains if any fault occurs. Pull out the plug, switch off or remove the fuse.
- 5) When disconnecting the appliance, pull on the plug, not on the cable.
- 6) Any repairs and work on the appliance should only be carried out by the customer service department, as unauthorised work could prove highly dangerous for the user. The same applies to changing the mains power cable.
- 7) Do not stand on the plinth, drawers or doors or use them to support anything else.
- 8) If you have a lockable appliance, do not keep the key near the appliance or within reach of children.
- 9) Do not store explosives or sprays using combustible propellants such as butane, propane, pentane etc. in the appliance. Electrical components might cause leaking gas to ignite. You may identify such sprays by the printed contents or a flame symbol.
- 10) The noise level while the appliance is operating is below 70db(A) (relative noise level 1 pW).

2. SETTING UP

- 1) Avoid positioning the appliance in direct sunlight or near cookers, radiators and similar sources of heat.
- 2) The floor on which the appliance stands should be horizontal and level. To compensate for any unevenness, adjust the height of the feet. To alter the height, rotate the bottom section of the respective adjustable foot

WARNING!


The height of the foot can be adjusted from 120 mm to 170 mm. Do not set the adjustable foot to a height greater than 170 mm! The bottom section of the adjustable foot can become loose and the appliance may then tip over. This can lead to serious or even fatal injuries.

- 3) There must be a gap of at least 30 cm between the upper edge of the appliance and the ceiling.

Do not cover ventilation openings or grille.

3. DIMENSIONS

Gross capacity (see type plate).

	Height	Width	Depth
TEMPLOW-L	2150-2200	700	830

4. ELECTRICAL CONNECTION

Only operate the appliance with alternating current (AC).


The permissible voltage and frequency are indicated on the type plate.

The socket must be properly earthed and protected by a fuse. The tripping current of the fuse must be between 10 A and 16 A.

The socket must not be situated behind the appliance and must be easily accessible.

Do not connect the appliance using an extension cable or extension socket.

Do not use stand-alone inverters (conversion of direct current to alternating/ three-phase current) or energy-saving plugs. Risk of damage to the electronic control system!

WARNING!

It is recommended to clean the appliance before switching it on for the first time (see "Cleaning" section).

Before switching the appliance on for the first time, plug the cleaning water drain opening in the floor of the appliance with the plug provided. This is important to ensure that the appliance works properly.


5. EQUIPMENT

Cleaning water drain opening

A drain hose with a R 3/4 connection can be fitted to the underside of the appliance.

The water which collects in the interior during cleaning can be drained off in this way.

An angled connector is supplied with the appliance.


Intermediate shelves

Place the supplied intermediate shelves onto the shelf rails.

WARNING

The maximum load per intermediate shelf is 20 kg.


Safety lock

The lock is equipped with a safety mechanism.

Locking the appliance

- Insert the key as shown by arrow 1.
- Turn the key 180° (2).


To unlock the appliance, the same procedure must be repeated in the same order.

6. CONTROL PANEL

-  On/Off button (switching the appliance on and off)
-  Reading the temperature progression
-  Selection buttons
-  SuperFrost button
-  Audible alarm Off button
-  Defrost button (for manually activating the defrost function)
-  Enter button

Symbols in the display

-  Compressor is running
-  LED flashing - refrigeration unit switches on after a delay. The compressor will start automatically after the pressure in the refrigerant circuit has equalised.
-  Fan is running
-  Appliance is defrosting
-  The  display means that the power supply and interior temperature of the appliance are recorded.
-  If  flashes in the display, there has either been a power failure or the temperature in the appliance exceeded the permissible range.
-  SuperFrost is activated
-  Alarm function
-  The appliance has suffered a fault. Contact the customer service department.

7. SWITCHING ON/OFF

Switching the appliance on

Press  for approx. 3 seconds. Display = **ON**

Switching the appliance off

Press  for approx. 3 seconds. Display = **OFF**

8. TEMPERATURE SETTINGS

Press  for 1 second. The temperature display flashes.

To increase the temperature (warmer): press button .

To reduce the temperature (colder): press button .

Press  again. The desired temperature setting is saved.


Note

The temperature in the warmest area of the interior may be higher than the temperature setting.

If the door is left open for a lengthy period, the temperature in the appliance's compartments may rise dramatically.

Temperatur display mode


The temperature display can be switched between degrees Celsius and degrees Fahrenheit. Factory setting is degrees Celsius.

Press  for 5 seconds. Display = 

Press . Display = 

Use button  or  to select the desired setting.

0 = °C 1 = °F



Press . Display = 

Press  for 5 seconds. The electronic control system will switch back to normal operating mode.

SuperFrost

Use SuperFrost, to rapidly cool large amounts of product. When SuperFrost is activated, the appliance operates with maximum cooling performance.

Activating SuperFrost

Press  for approx. 3 seconds. Display =   


Place the products inside the appliance.

The electronic control system will automatically switch back to normal operating mode.

9. ALARMS

Door open alarm

When the door is opened, the LED  lights up and the temperature display begins to flash.

When the door has been left open for more than 240 seconds, the LED  begins to flash, and  and the temperature indication flash alternately in the display.


The audible warning signal sounds (unless the audible warning signal function has been deactivated).

If the door has to stay open for longer in order to insert items to be cooled, cancel the audible warning signal by pressing button .


Setting the delay time for the door open alarm

The time before the audible warning signal sounds after the door has been opened can be adjusted.

Press  for 5 seconds. Display = r15

Press  until   appears in the display.

Press . Display = | Setting range = 1 - 5 minutes.

Use button  or  to select the desired setting.

Press . Display =  

Press  for 5 seconds. The electronic control system will switch back to normal operating mode.


Deactivating the audible warning signal function

The audible warning signal function can be completely deactivated if necessary.

Press  for 5 seconds. Display = r15

Press  until   appears in the display.

Press . Display = 

Use button  or  to select the desired setting.

0 = activated 1 = deactivated

Press . Display =  

Press  for 5 seconds. The electronic control system will switch back to normal operating mode.

Alarm test

This test checks the function of the internal and any external connected alarm device. The appliance does not stop its refrigerating function during this test.

Press  +  for 5 seconds.

- The display will change to a temperature value of 0.2°C below the set upper alarm limit.
- The temperature value will now rise by 0.1°C every 2 seconds.
- When the upper alarm limit is reached, **HIO** will appear in the display. An external alarm unit connected to the floating alarm output will now be activated.
- The temperature value will continue to rise up to 0.2°C above the upper alarm limit.
- The same process will take place automatically for the lower alarm limit. **LIO** will appear in the display.

The LED  will be lit during the test.

The electronic control system will automatically switch back to normal operating mode.

Cancelling the test prematurely

Press  for 5 seconds.

Note

If the values of the upper and lower alarm limit (**AL** and **AH** in the section entitled "Adjusting the alarm parameters") are set to **0**, **H--** and **L--** will appear in the display during this test.

Alarm messages

1. LED flashes in the display

If  appears in the display, the appliance has a fault. Consult your nearest customer service point.

2. LED flashes in the display; the display reads HI or LO

The interior is too warm (HI) or too cold (LO).

The audible warning signal sounds (unless the audible warning signal function has been deactivated).

Note

The alarm parameters can be adjusted. See **Adjusting the alarm parameters**.

3. HA / flashes in the display


The interior was too warm or too cold (**HA**) during a certain period of time.

Adjusting the alarm parameters


The alarm limits (difference to the set temperature) and the alarm delay (delay until alarm sounds) can be adjusted.

Press  for 5 seconds. Display = r15


Press  until   appears in the display.


  = Lower alarm limit

Press . Display = temperature difference in °C


Use button  or  to select the desired setting.

Set positive values only.


Press . Display =  


Press . Display =   Upper alarm limit

Press . Display = temperature difference in °C


Use button  or  to select the desired setting.


Set positive values only.


Press . Display =  

Press . Display =  

Press . Display = alarm delay in minutes

Use button  or  to select the desired setting.

Press . Display =  

Press  for 5 seconds.

The electronic control system will switch back to normal operating mode.

Reading the temperature progression

Press $\overline{\text{HACCP}}$. Display = $r \text{ } \bar{t}$

Scroll through the list using ∇ or \wedge .

$r \text{ } \bar{t}$ Period in hours in which the maximum and minimum interior temperatures were measured

$r \text{ } H$ Maximum (highest) measured temperature

$r \text{ } L$ Lowest measured temperature

Select the required item using the ⚙ button. Press this button again to return to the list.

You can exit the menu at any time by pressing ⏏ for 5 seconds.

If no button is pressed within 60 seconds, the electronic control system switches back automatically.

Resetting the alarm events

Press $\overline{\text{HACCP}}$. Display = $r \text{ } \bar{t}$

Press $\overline{\text{HACCP}}$ + \wedge for 5 seconds. Display = $r \text{ } \bar{E} \bar{S}$.

Press ⏏ for 5 seconds. The electronic control system will switch back to normal operating mode.

Resetting the recorded temperature progression $r \bar{t}$

Press $\overline{\text{HACCP}}$. Display = $r \text{ } \bar{t}$

Press ⚙ . Display = $\square \cdot 999$

Press ∇ for 5 seconds. Display = $r \text{ } \bar{E} \bar{S}$.

The values for $r \text{ } H$ and $r \text{ } L$ (highest and lowest measured interior temperature) are then reset to the current interior temperature.

Press ⏏ for 5 seconds. The electronic control system will switch back to normal operating mode.

Resetting the parameters to factory settings

The alarm limits and sensor calibration values can be reset to the factory settings using this function.

Pull out the mains plug.

Keep ⏏ pressed and connect the mains plug.

Display = $\bar{b} \bar{r} \bar{1}$

Press ⚙ . Display = $\bar{5} \bar{t} \bar{d}$

The electronic control system will switch back to normal operating mode.


Storing

The grid shelves can be moved to accommodate different height bottles or packages.

(1) The ventilation slots on the inside must not be covered!

(2) Stacking mark ⏏

Only load the top shelf up to the stacking mark. This is important so as to ensure that the air can circulate properly and the temperature is even throughout the interior.

NOTICE

Use the bottom shelf of the refrigerator/freezer compartment so that liquids cannot drip onto other products.

10. DEFROSTING

The refrigerator compartment defrosts automatically.

Activating the defrost function manually

If the door has been left slightly open for a long time, a layer of ice may form in the interior and on the cooling plate. The defrost function can then be activated manually.

Press ❄ for 3 seconds. Display = $\text{❄} + \bar{d} \bar{F} \bar{b}$

The electronic control system will automatically switch back to normal operating mode.

Display = $\bar{d} \bar{F} \bar{E}$

Setting the display indication for the defrost phase

Press ⏏ for 5 seconds. Display = $r \text{ } \bar{1} \bar{5}$

Press \wedge until $\bar{d} \bar{F}$ appears in the display.

Press ⚙ . Display = $\bar{1}$

Use button ∇ or \wedge to select the desired setting.

0 = Symbol ❄ + alternating display of $\bar{d} \bar{F} \bar{F}$ and the current temperature in the interior of the appliance.

1 = Symbol ❄ + temperature before the start of the defrost phase.

2 = Symbol ❄ + $\bar{d} \bar{F} \bar{F}$.

Press ⚙ . Display = $\bar{d} \bar{F} \bar{E}$

Press ⏏ for 5 seconds. The electronic control system will switch back to normal operating mode.

11. CLEANING

CLEAN THE APPLIANCE AT LEAST TWICE PER YEAR!

Before cleaning always switch off the appliance. Pull out the mains plug or switch off or unscrew the fuse.


- Empty the appliance and store the product in a cool place.
- Clean the inside and equipment with lukewarm water and a little detergent. Do not use abrasive or acid cleaners or chemical solvents.

Do not use steam cleaners because of the risk of injury and damage.

- Ensure that no cleaning water penetrates into the electrical components or ventilation grille.
- Dry all parts well with a cloth.
- Use a commercially available stainless-steel cleaning agent for stainless-steel appliances.

Do not use abrasive sponges or scourers, do not use concentrated cleaning agents and never use cleaning agents containing sand, chloride or acid or chemical solvents, as these would damage the surfaces and could cause corrosion.

Neutral PH multi-purpose cleaners are recommended.

Cleaning the dust filter

Clean the dust filter at least twice per year!
Remove the plug!

1. Remove the dust filter by lifting upwards.
2. Clean the dust filter with water and detergent.
3. Reinstall the dust filter.

12. SHUTTING YOUR APPLIANCE DOWN

If your appliance is to be shut down for any length of time, switch it off and disconnect the plug or switch off or unscrew the fuse.

Clean the appliance and leave the door open in order to prevent unpleasant smells.

The appliance complies with the relevant safety regulations and EC Directives 2004/108/EC and 2006/95/EC.

13. MALFUNCTIONS

You may be able to rectify the following faults by checking the possible causes yourself:

- Appliance does not function:
 - Is the appliance switched on?
 - Is the plug correctly fitted in the mains socket?
 - Is the fuse intact?
- The temperature is not low enough:
 - Is the temperature setting correct (see "Setting the temperature")?
 - Have excessive amount of products been placed in the appliance?
 - Does the separately installed thermometer show the correct reading?
 - Is the ventilation system working properly?
 - Is the appliance set up too close to a heat source?


If none of the above causes apply and you cannot rectify the fault yourself, contact the technical service department stating the type designation and serial number as indicated on the type plate.

Possible error messages in the display

Error code	Error	Action
E0, E1, E2, rE	Temperature sensor defective	Contact the customer service department
EE, EF	Electronic control system error	Contact the customer service department
dOr	Appliance door open for too long	Close appliance door
HI	Temperature inside appliance too high (too warm)	Check that the door has been closed properly. If the temperature does not drop, contact the customer service department.
LO	Temperature inside appliance too low (too cold)	Contact the customer service department
HA	The interior was too warm or too cold during a certain period of time.	See Reading the temperature progression

14. CHANGING OVER DOOR HINGES

Door hinges should only be changed by a trained expert.
Changing the door hinges must be done by two people.

1. Open door by about 90°.


NOTE

The door must be opened 90° before the lower hinge bracket is removed.

This will hold the self-closing mechanism that is integrated into the door in the required position for installation.

WARNING!


If the door is removed and reinstalled in the closed position, this will lead to destruction of the self-closing mechanism on the first opening of the door.


2. Remove hinge bracket. Pull the door out at the bottom and lift off.

IMPORTANT


As a result of its heavy weight, the door must be held secure by one person.


3. Unscrew base mounting bracket.

4. Re-fit base mounting bracket on the opposite side.


5. Remove bearing ring from upper door mounting.

6. Remove the hinge components from the hinge bracket.


Note

The hinge bushing with brass washer and bearing ring may stick in the door mounting when you pull out the hinge bracket and, in this case, must be removed from it.


7. Transfer the upper hinge bracket and covers to the opposite side.

8. Turn the door by 180°.

9. Insert the bearing ring and hinge bushing in the upper door mounting.

10. Keeping door open at 90°, suspend in top square pin.


WARNING!

It is essential that the door is open at an angle of 90° during installation.


If the door is installed in the closed position, this will lead to destruction of the self-closing mechanism on the first opening and closing of the door.

IMPORTANT


As a result of its heavy weight, the door must be held secure by one person.

11. Fit the hinge components on the hinge bracket.


12. Install the hinge bracket with the fitted hinge components in the lower door mounting.

Note

The tabs on the bearing ring must fit into the recesses of the door mounting during installation.

If necessary, gently rotate the upper ring.

13. Screw hinge bracket into place.