**pH-metro digital portatil**

**Portable digital pH-meter**

**pH-2004**

**4120400**


# Indice

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# Table of contents

<b>1 Packing list and accessories</b> .....	<b>12</b>
1.1 Packing list .....	12
1.2 Accesories .....	12
<b>2 Safety</b> .....	<b>13</b>
<b>3 Overview</b> .....	<b>13</b>
<b>4 Technical specification</b> .....	<b>13</b>
<b>5 Description</b> .....	<b>14</b>
5.1 Display .....	14
5.2 Rear panel .....	14
<b>6 Operation</b> .....	<b>14</b>
6.1 Placing batteries .....	14
6.2 pH electrode and temperature probe connection .....	14
6.3 Panel .....	14
6.4 Initial screen. Measuring pH or mV. ....	16
6.5 Printing the readings .....	16
6.6 Storing readings .....	16
6.7 Deleting readings .....	16
6.8 Viewing stored readings .....	16
6.9 Calibrating pH-2004 .....	17
6.10 Isopotencial point selection .....	18
6.11 Measuring temperature. Automatic compensation. ....	18
6.12 Automatic switching Off .....	18
6.13 Buffer solutions. ....	18
6.14 RS-232 Communication .....	18
6.15 Electrode handling .....	19
<b>7 Warranty</b> .....	<b>19</b>
<b>8 Statement of conformity CE</b> .....	<b>19</b>

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

  
David Pecamins  
Rble. Calidad

## 1 Packing list and accessories

First, check that the content is complete and looks in a good state.  
Keep the transportation box for some days.

### 1.1 Packing list

- pH-meter (1)
- Arm with probe holder (2)
- AA size, 1.5 Volt batteries (4 units) (3)
- Short-circuited connector. (4)
- RS-232 cable (9) (5)
- RS-232 cable (20) (6)
- Case (8)
- pH electrode + temperature probe (9)
- Instruction manual

### 1.2 Accesories

#### pH + Temperature probe

Code4120401

Set of pH electrode (Epoxy) and metal Rt 2252 temperature probe. Connectors are IP65 sealed.

#### Probe holder

Code1001552

Up to 4 probe on the same holder.  
Multiposition arm.

#### Magnetic stirrer «agimicro»

Code 5504300

Regulable speed from 0 to 1400 r.p.m.  
Engineering material chemical resistant.  
Included, magnet PTFE coated Ø8 x 32 mm.

#### Certified Buffer solutions for calibration

Buffer pH=4.00 (20°C)    pH=4.01 (25°C)  
Tampón pH=6.88 (20°C)    pH=6.86 (25°C)  
Buffer pH=9.23 (20°C)    pH=9.18 (25°C)  
Electrolyte KCl 3M

Code 4120107  
Code Please ask.  
Code 4120120  
Code 4120109

#### Serial printer

Code 4120113

Serial Ink ribbon printer. (Not thermal paper)  
Direct connection to pH-2006.  
Paper roll 2 1/4"

Warning symbol

Important  
Information  
symbol

## 2 Safety

Across this manual, some risk of danger situations are remarked. Proceed as described to avoid any hazard.

- **Warning symbol of a dangerous operation:** Proceed according indicated instructions.
- **Important Information symbol:** Information important to obtain good results or for an optimal operation of the equipment or to extend the life of the equipment or to avoid the degradation of some of its components.
- Using a pH-meter, the operator, could be exposed to chemical hazard.
- This equipment must only operated by personnel trained an suitable qualified to work in a chemical laboratory and its common hazards and who have read and understood this manual or have been qualified on its manipulation and operation.
- This equipment its not designed and not manufactured to work in an explosive ambient area.

## 3 Overview.


The pH-2004 is a portable pH-meter protected with a water proof degree. Suitable for measurements outdoors (lakes, rivers,..) or on the laboratory.

The pH-2004 also measures temperature for ATC. If no temperature probe is connected, sample temperature can be entered by using the keypad and also ATC is applied on the measure. (ATC: Automatic temperature compensation)

The reading in can be used to take measures with Ion selective and Redox probes.

## 4 Technical specification

Power supply:	4x 1.5V AA size. Alkaline batteries
pH range:	0..14
pH precision:	± 0.01 pH
pH resolution:	0.01 pH
mV range:	±1800 mV
mV precision:	±1 mV
mV resolution:	1 mV
Temperature correction range:	0 ... 60 °C
Temperature precision:	±0.5°C
Temperature resolution:	0.1°C
Temperature compensation precision:	±0.015 pH
Input impedance:	1 x 10 <sup>12</sup> Ohm
Calibration points:	2
Ambient working conditions:	5 ... 40°C humidity (Hr) lower than 80%
Dimensions:	201 x 100 x 45 mm
Weigth:	0.5Kg
Protection degree against water:	IP 65
Temperature probe type:	Rt-2252

Electrode pH    Temperature probe    Output RS-232

pH electrode and temperature probe.

Electrode cover.

## 5 Description

All function units are selected from the front keypad. The display shows the pH reading, the temperature and the pH-meter status.

### 5.1 Display

- pH or mV reading.
- Units (pH / mV).
- Measuring temperature.
- Status indicator.

### 5.2 Rear panel

- BNC socket for pH electrode.
- Mini BNC socket for the temperature probe.
- Serial RS-232 socket.

## 6 Operation

Once unit is unpacked check packing list. Visually check that all goods. Keep the package for some days.

### 6.1 Placing batteries

- Place the batteries on the case. (Rear side)
- Pay attention on polarity.
- Replace batteries when «LoBat» appears on the display.

### 6.2 pH electrode and temperature probe connection

The electrode must be stored with its immersion cap filled with electrolyte. The electrode and the temperature probe must be connected on the rear panel. See 5.2.

To connect the electrode, remove the short-circuit connector from the pH-meter and to connect the electrode to BNC socket. Please keep the BNC socket with the short-circuit connector when the pH electrode is not connected on the socket.

**Note:** The BNC connectors are extracted and placed turning 1/4 turn.

Keep the electrode with its immersion cap until use.

### CAUTION:

The immersion cap contend KCl 3M as electrolyte which is corrosive.

Take suitable cautions to avoid electrolyte be in contact with skin or wear. Avoid any contact with eyes.

### 6.3 Panel

Switch on the pH-meter by pressing «ON/OFF» key.

Wait until the unit act before to press any other key.

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




Press «ENTER» to view next reading.

#### 6.4 Initial screen. Measuring pH or mV.

After switch on pH-2004 (Key «ON/OFF») and probes are connected, it starts to read pH and temperature. Press «▲ /pH/mV» to read mV.

The reading in mV can be used to work with ion selective electrode or REDOX electrode. It is also useful for electrode status evaluation.

#### 6.5 Printing the readings

If the pH-2004 has been connected to a serial printer(RS-232) reading can be printed.

- Instantaneous printing: Pulsar «ENTER / PRINT»
- Printing stored readings: (See 6.6).
- First select pH or mV. Because both are printed separately.
- Press «MODE/MEAS»
- Press «▼ / SAVE» until get PRN screen.
- Press «ENTER/PRINT» to validate.

#### 6.6 Storing readings

Directly from the measuring screen press «▼ / SAVE» to store the actual reading. The pH and mV reading are stored, printed and deleted separately.

First select pH or mV.

Up to 250 reading can be stored. If more than 250 are stored the new ones are stored over the oldest.

#### 6.7 Deleting readings

Deleting stored reading:

- First select pH or mV. Because both are stored separately.
- Press «MODE/MEAS»
- Press «▼ / SAVE» to DEL screen.
- Press «ENTER/PRINT» to validate.

#### 6.8 Viewing stored readings

Viewing stored readings:

- First select pH or mV. Because both are viewed separately.
- Press «MODE/MEAS»
- Press «▼ / SAVE» to VIEW screen.
- Press «ENTER/PRINT» to validate.

### 6.9 Calibrating pH-2004


The pH-2004 need the following standards for calibration: (See 1.2)

Starting screen for calibration.

T <sup>a</sup>	pH 1	pH 4	pH 7	pH 9	pH 12
20°C	1.675	4.002	6.881	9.225	12.627
25°C	1.679	4.008	6.865	9.180	12.454
30°C	1.683	4.015	6.853	9.139	12.289

Calibrate in one or two points using one or two of this buffers.

Starting screen for one point calibration or first point calibration.

Choose this two point by considering:

- One point is always pH=6.85
- If you are going to measure an acid solution choose pH=4.00
- If you are going to measure an alkali solution choose pH=9.18

Starting screen or second point calibration.

Calibrating pH-2004:

- Clean the electrode with water.
- From the measuring screen, press «MODE/MEAS»
- Press «▼» to get «STD» screen validate with «ENTER».
- Immerse the electrode on the buffer solution pH=6.86 at 25°C
- Wait for stable pH reading.  
Reading can be switched to mV by key "pH/mV".
- Press «ENTER»
- The display shows the recognized buffer value.  
This value can not be changed.  
After 5s. The unit jumps, automatically to STD2 screen.
- Press «ENTER»

If calibration its only in 1 point, that's all, press «MODE/MEAS»

To continue for a two points calibration (strongly suggested) press «ENTER» on «STD2» status.

- Clean, rinse, the electrode with distilled water.
- Immerse the electrode on the buffer solution pH=4.00 or 9.18 a 25°C
- Wait for stable pH reading.  
Reading can be switched to mV by key "pH/mV".
- Press «ENTER»
- The display shows the recognized buffer value.  
This value can not be changed.  
After 5s. ...
- The calculated SLOPE value is shown

This value can be used to evaluate the electrode state:

Slope:	electrode state:
98-105%	Electrode is good.
94-97%	Clean the electrode with water and calibrate again.
90-93%	Bat state /repair or change the electrode. See 6.15


**Nota:** Calibrate pH-2004 when:

- First time use
- Daily.
- After change the electrode
- After measuring strong acids (pH<2)
- After measuring strong alkalis (pH>12)
- After measure very dense solutions that can obstruct the membrane.



Isopotential point selection screen

Entering measuring temperature screen.

Enter measuring temperature.

Automatic pH-meter disconnection.

### 6.10 Isopotential point selection

Select suitable isopotential point to measure pH value of pure water, super pure water and boiler water.

Usually, isopotential point is set to 7.00pH to measure general water solution.

To measure pure or super pure water set 12.00pH

For ammonia water set 17.00pH.

### 6.11 Measuring temperature. Automatic compensation.

The pH-2004 includes ATC (Automatic temperature correction). This means that the reading is corrected if the measuring temperature is different from 25°C.

The measuring temperature is the one the temperature probe reads or if the temperature probe is NOT connected it is the temperature set by the keypad.

Setting the measuring temperature by the keypad:  
(Temperature probe NOT connected on the rear panel socket)

- Press «MODE/MEAS» and select «MTC» using arrow keys.
- Validate with «ENTER».
- Use keys «▼» «▲» to select the temperature.
- Validate with «ENTER» and press «MODE/MEAS» to return.

**Note:** If the measuring temperature set have a value far from the actual sample temperature the pH reading could have considerable error.

### 6.12 Automatic switching Off

Sets the automatic disconnection after 2h of no operation.

- Press «ENTER».
- Use keys «▼»: y «▲»  
ON : Automatic switch Off  
OFF: No automatic switch Off

### 6.13 Buffer solutions.

For a certified calibration, some standard buffers must be used, As the ones described on 1.2.

Clean, rinse, with water, the electrode before to immerse into buffers.

#### Attention

- Buffer solutions are irritant chemicals. Avoid any contact with eyes and use a lab coat to protect your wear.
- In case of chemicals be in contact with skin, clean, quickly, with abundant water.
- To prepare the buffer solutions just take a 300-400 ml flask, and solve the poulder bag into 250ml of distilled water. Stirring until solved.
- Keep this solutions in a closed flask. Do not use this buffer for calibration after 1 moth after preparation.

### 6.14 RS-232 Communication

Use the communication software to capture reading on the computer. Connect the communication wire on a free RS-232 COM port on your computer. (Speed 9600)

### 6.15 Electrode handling

- If the electrode comes with his own instruction sheet, follow it.
- Keep the electrode with its end cap filled with KCl 3M when not used.
- Take suitable cautions to avoid electrolyte KCl 3M be in contact with skin or wear. Avoid any contact with eyes.
- Avoid than the electrode left more than 1h without the end cap and becomes dry on the air. It can lost its features if becomes dry. For a short periods left it immersed on a glass of water.
- Clean and rinse with abundant water after measure.

## 7 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep purchasing bill as the only valid document for the purchasing date.

The warranty does not covers the damages caused for a unit malfunction. The unit must be suitably qualified before to be used for valuable samples processing.

Return good are not accepted without previous J.P. SELECTA, s.a. approval.

## 8 Statement of conformity CE

The manufacturer:

**J.P. SELECTA, s.a.**  
Ctra. NII Km 585.1  
08760 ABRERA (BARCELONA)  
ESPAÑA

Declares that the laboratory equipment:

pH-meter pH-2004 Code: 4120400

Meet the following European Directives:

73/23/CEE Directive for electrical safety.  
89/336/CEE Directive for electromagnetic compatibility.(EMC)

Meet the following Standards:

EN 61326:2002  
Electrical equipment for measurement, control and laboratory use. EMC Requirements.

EN 61010-1:2002  
Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1: General requirements.

Abrera, June 2006

Ramón Ramón  
Dtor. Técnico.

David Pecanins  
Rble. Calidad

J.P. SELECTA, s.a.

---