HOTCOLD F
HOTCOLD J

2101517
2101518

ESTUFA REFRIGERADA
REFRIGERATED CABINET

INDICE

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GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts, please inquiries contacting with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 2009/104/CE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

PACKING LIST

The standard equipment consists of the following components:

Description	Code	Quantity
Hotcold F	2101517	1
Hotcold J	2101518	1
Shelves	1001858	3
Instruction Manual	80411	1

ACCESSORIES

Hotcold F and J equipments allow to connect stirrers and shakers without heating elements, which can be connected on inner sockets or through lateral hole.

TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Code		2101517	2101518
Temperature range		+6°C to 20°C	+6°C to 20°C
Humidity range		26% to 33%	47% to 60%
Capacity (litres)		600	600
Inner useful dimensions (cm) *with / without shelf	Height	138	138
	Width	53/61*	53/61*
	Depth	58/68*	58/68*
Overall dimensions (cm)	Height	210	210
	Width	73	73
	Depth	89	89
Shelves (maximum)		18	18
Power (HP)		3/4	3/4
Consumption (W)		1350	1350
Weight (kg)		133	133

EQUIPMENT DESCRIPTION

Hotcold F and J refrigerated cabinets are built in AISI 304 stainless steel.

These are fan air circulation equipments which allow:

Stability*:	±0.5°C
Homogeneity*:	±1.3°C
Resolution:	±0.1°C
Humidity:	±6%

*Note: tests carried out in environmental conditions 20°C and 50% RH

These equipments are designed for working under laboratory conditions of:

Temperature	18-27°C
Relative humidity	25-75%

Working outside these margins can limit their benefits.

Note: Humidity is a fixed value on settings. It is not possible for the user to modify it.

These cabinets are commonly used to preserve products under fixed conditions. The equipments do not generate humidity. Very dry environments can affect the operation of these equipments.

INSTALLATION

Place the cabinet on a flat horizontal, level surface, leaving a space of about 10 cm at the back and sides of the equipment.

For a maximum efficiency, do not place the machine in direct sunlight or other sources of heat.

Connect the equipment to the mains. Be sure that the voltage supply is the same as the one indicated on the characteristics plate. Do not use the equipment if it is NOT earthed.

Note: Before connecting it to the mains, it is advisable to leave the machine at rest for 24h to stabilize oil and gas cooling unit.

If you have to change the plug cap, bear in mind the following:

Blue cable: Neutral.

Brown cable: Phase.


Yellow/green cable: Earth.

Check that the circuit breaker on the electrical panel (4) is powered. See "Electrical panel" section on the "Operation" chapter.

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
Do not use the cabinet for thermal treatment of products which can give off vapours capable of making an explosive mixture.

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

Inside view of the cabinet

INNER CHAMBER:

- 1. Internal light.
- 2. Air circulating fan
- 3. Side hole for connections inlet.

**ELECTRICAL PANEL:
(back upper side of the equipment)**

- 4. Magnetothermic control, cold - heating group and plugs.
- 5. Safety thermostat


OPERATION

All equipment functions are set using the touch screen. Gently press your finger in the center of each icon.

1 – CONTROL PANEL

The control panel is managed by using virtual keys. Press gently on the key that appears on the screen, as closed as possible to its center.

Setpoint temperature
Inside chamber temperature

Relative humidity value indicator
Technical service menu key

To modify the setpoint temperature value, click on the setpoint value. It will change color. Press on the increase/decrease keys to adjust the temperature value.

Once the value is adjusted, the system will automatically change the setpoint value.

2 – TECHNICAL SERVICE MENU

Press the enter key and the system will ask for a password to enter the setup screen. EXCLUSIVE USE FOR TECHNICAL SERVICE.


To enter this menu, the system will ask for a PIN code (0284)

Click on each X and modify its value with the arrow keys, to enter the correct pin.

Press the return key to exit this screen. Press the OK key to access the next screen.

If the pin is correct, the following screen will appear:

1. Time between defrosts. Factory value: 8 hours.
2. Defrost time. Factory value: 5 minutes.
3. Fan start timeout. After a defrost cycle, to prevent heating from moving inside the chamber, the fan waits this time. Factory value: 1min
4. Temperature offset correction. Adjust this value by comparing with a standard thermometer.
5. Maximum defrost temperature. The internal temperature probe measures the temperature of the defrost resistance, if this value is exceeded, it cuts off the power supply to this resistance. Factory value: 80°C.
6. Relative Humidity offset correction. Adjust this value by comparing with a standard hygrometer.

OK key to accept the values.

RETURN key to exit without saving the values.

STARTING UP:

1. Check that the installation protocol has been correctly followed.
2. Switch on the main switch (see picture) located on the back side of the equipment. The graphical display of the temperature controller will light up.

MAINTENANCE

Before removing the case, disconnect the equipment from the mains.

The manipulation of the internal electronic circuits of the cabinet by unauthorized personnel can cause irreparable damages. Take it to one of J.P.SELECTA's authorized technical services.

CLEANING:

For cleaning different parts, we recommend to use the following products:

Stainless steel cleaning: Alcohol

Plastic cleaning: Alcohol with cotton duster or with a non-abrasive cloth.

SPARE PARTS

16099	Temperature regulator
20075	Magnetothermal
29553	Control Circuit
36043	Programmable relay
36027	Solid relay 25A
39057	Electrical resistance
43031	Safety thermostat
43071	Temperature probe
85167	4.3" TFT screen


GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.

SAFETY THERMOSTAT

- 8. Safety thermostat regulation screw (grey).
- 9. Safety thermostat reset button (white).

Rearming the safety thermostat

Safety thermostat is an element to prevent the temperature from exceeding a value. If by a fortuitous event the thermostat is triggered, that is, it opens the power circuit, the heating element will stop working. Usually there is an indicator light, or an error message, that indicates that the thermostat is in this position.

If the thermostat is triggered, check that the equipment looks good and it is not overheating at any point.

To change its status, simply press the white button (9), and you will hear a "click" that indicates that the thermostat has been reset. If the equipment has not cooled down sufficiently, the thermostat cannot be reset, and you must wait for the temperature to drop before it can be reset.

Adjusting the safety thermostat

From factory, the thermostat is adjusted to its maximum temperature. To adjust the safety thermostat temperature, do as follows:

1. With the equipment connected to the mains, turn the safety thermostat screw (8) clockwise to the end of the scale, with the help of a flat screwdriver.
2. Select the desired safety temperature using the control keys. This temperature has to be minimum at 3°C above the working temperature.
3. Turn on the equipment with the Start key.
4. When the equipment is stabilized at the above selected temperature, turn the safety thermostat screw (8) slowly to the left until it clicks, and the red indicator lamp lights up.
5. Stop the equipment and leave it to cool with the lid open for 15 minutes.
6. Select the working temperature by using the temperature selection keys (minimum at 3°C below the desired safety temperature).
7. Turn the equipment on again with the Start key. Once the temperature is stabilized, the machine is ready to work.