PHARMALOW - M

2101287

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1. SAFETY INSTRUCTIONS

- 1) To prevent injury or damage to the unit, the appliance should be unpacked and set up by two people.
- 2) In the event that the appliance is damaged on delivery, contact the supplier immediately before connecting to the mains.
- 3) To guarantee safe operation, ensure that the appliance is set up and connected as described in these operating instructions.
- 4) Disconnect the appliance from the mains if any fail occurs. Pull out the plug, switch off or remove the fuse.
- 5) When disconnecting the appliance, pull on the plug, not on the cable.
- 6) **Any repairs and work on the appliance should only be carried out by the customer service department, as unauthorized work could prove highly dangerous for the user.** The same applies to changing the mains power cable.
- 7) Do not stand on the plinth, drawers or doors or use them to support anything else.
- 8) **If you have a lockable appliance, do not keep the key near the appliance or within reach of children.**
- 9) Do not store explosives or sprays using combustible propellants such as butane, propane, pentane etc. in the appliance. Electrical components might cause leaking gas to ignite. You may identify such sprays by the printed contents or a flame symbol.

2. SETTING UP

- 1) **Avoid positioning the appliance in direct sunlight or near cookers, radiators and similar sources of heat.**
- 2) The floor on which the appliance stands should be horizontal and level. Compensate for uneven floors with the adjustable feet. Leave a minimum space of 10cm between the back of the equipment and the wall.
- 3) Always ensure that there is a good ventilation and that the outward flowing air is able to escape.

3. FEATURES

- Gross capacity (see plate): 360L
- External dimensions: 185 x 60 x 62cm (H x W x D)
- Internal dimensions: 145 x 47 x 55cm (H x W x D)
- Homogeneity: 1.5°C
- Stability: 1°C
- Resolution 1°C

4. CONDITIONS OF USE

- Altitude up to 2000 m (795 mbar)
- Temperature from 12°C to 32°C (53,6°F to 89,6°F)
- Maximum relative humidity 65 %
- Mains voltage fluctuations up to $\pm 10\%$ of rated voltage.
- Equipment for indoor use.

5. POWER CONNECTION

Power supply (a.c.) and voltage at the operating point must comply with the details on the type plate, which is attached at the rear part of the refrigerator. The socket must be fused with a 10 A fuse or higher.

To ensure trouble-free operation, do not connect the appliance to the power supply until 6 to 8 hours after delivery. Connect the appliance through a properly earthed plug and socket.

The socket should be easily accessible and not directly behind the appliance.

The wires in the mains lead are coloured in accordance with the following code: green/yellow = earth, blue = neutral, brown = live.

Warning! This appliance must be earthed.

6. SWITCHING THE APPLIANCE ON AND OFF

You are recommended to clean the appliance before switching it on for the first time (see "Cleaning"). When it is connected up (plugged into an approved wall socket) the appliance is ready for operation.

To switch the appliance off: pull out the mains plug.

7. OPERATING

Our products leave the manufacturing process after through testing procedures to ensure they are in proper working order.

To operate them, simply connect them to a GROUNDED electrical outlet to the proper voltage of the equipment. The unit should operate normally.

The tubes are inside the chamber. This switch should not affect any other components.

An electronic alarm is also incorporated and it advises when the chamber temperature reaches dangerous levels (up or down) for correct conservation of the medicines.

The alarm is visual and shows "LA" for temperatures under minimum or "HA" for programmed maximum temperatures.

Immediately, make sure that the COMPRESSOR and the FAN begin to function normally (they are located at the rear part of the unit).

Also check that the TURBINE on the upper part of the chamber spins normally. From this point forward the unit functions automatically; the temperature is controlled by the electronic THERMOSTAT THERMOMETER which regulates on and off in order to maintain the unit within the required temperature range (stops at +3°C and starts at +6°C). If no malfunction is noted to this point, allow the unit to function EMPTY during TWO or THREE HOURS AT LEAST so it can reach NORMAL OPERATING CONDITIONS.

We have programmed an automatic defrosting cycle approximately every four hours to eliminate ice build-up on the evaporator. This frost could form an icy film that could damage the evaporator if it is not removed in time.

This defrosting cycle lasts approximately TEN MINUTES. During this time, the compressor stops working, interrupting the refrigeration cycle in order to remove the ice build-up on the evaporator.

During this cycle we will notice a rise in the temperature inside the chamber; however, this rise in temperature is normal: Given the natural THERMAL INERTIA, a longer stoppage time and steeper rise in air temperature would be required in order to negatively affect its INTERNAL TEMPERATURE.

At the end of the defrosting cycle, the temperature recovers rapidly, without causing damage to the articles in storage.

The defrosting water is automatically evaporated. No manual defrosting is required.

Do not store articles in the refrigeration unit until all of the instructions noted above have been completed. The unit must be in normal operating conditions before items can be stored in it.


We recommend periodic observation of the condenser unit. If dirt build-up is noticed, it should be cleaned with a brush.

We strongly recommend against manipulation of the components of this unit by non-authorized personnel.

If you observe any malfunction, contact technical service immediately.

8. TEMPERATURE SETTINGS

The following icons appear on the screen:

LED	ON	OFF	FLASHING
	compressor ON	compressor OFF	- Compressor protection activated - Setpoint setting
	Evaporator fan ON	Evaporator fan OFF	Fan shutdown activated
	Interior light ON	Interior light OFF	Interior light activated via digital input
AUX 1	Auxiliary function 1 ON	Auxiliary function 1 OFF	- Auxiliary function 1 via digital input - Auxiliary function 1 via digital input
AUX 2	Auxiliary function 2 ON	Auxiliary function 2 OFF	- Auxiliary function 2 via digital input - Auxiliary function 2 - delay activated
	Defrost or pre-dripping activated	-	- Defrost delay activated - Dripping time enabled
	- Energy saving activated - Low power activated	-	-
	View time	-	date, time and day of the current week setting
	View temperature	-	Fast refrigeration activated
HACCP	recording of HACCP alarms	-	new HACCP alarm saved
	Alarm activated	-	-

Press any key for 1 second to unlock the keyboard and the screen shows the unlocking message 'UnL'.

Some functions that can be performed

To change the Set Point, check first that the keyboard is unlocked.

	Press the SET key.
	Press the up or down key to set a value within the r1 and r2 limits.
	Press the Set key (or do not press for 15 s).

For a manual defrost, check that the keyboard is unlocked, and that quick refrigeration is not active.

	Press the Defrost key for 2 s.
---	--------------------------------

If P3 = 1 (by default), defrost is activated whenever the evaporator temperature is below limit d2.

Changing parameters:

	Press the SET key for 4 seconds: the screen will show the message 'PA'.
	Press the SET key.
	Press the up and down keys to enter the password (by default '-19').
	Press the Set Key (or do not press for 15 s). The screen will display the message 'SP'.
	Press the up and down keys to select a parameter.
	Press the SET key.
	Press the up and down keys to select a value.
	Press the Set key (or do not press for 15 s).
	Press the SET button for 4 seconds (or do not press for 60 s) to end the process.

If the device is turned on, the screen will display the value P5 ('Temperature inside chamber'); if the screen displays any alarm code, go to section '**Alarms**'.

	If POF = 1 (by default), press the ON/STAND-BY key for 2 seconds.
---	---

Alarms

CODE.	DESCRIPTION	RESTART	TO CORRECT
Pr1	Chamber probe alarm	automatic	- check P0
Pr2	Evaporator probe alarm	automatic	- check the status of the probe
Pr3	Auxiliary probe alarm	automatic	- check the electrical connection
rtc	Clock alarm	manual	Set date, time and day of the week
AL	Low temperature alarm	automatic	check A0, A1 and A2
AH	High temperature alarm	automatic	check A4 and A5
id	Door open alarm	automatic	check i0 and i1
PF	Power fault alarm	manual	- Press a key - Check the electrical connection
COH	High condensation warning	automatic	check C6
CSd	High condensation alarm	manual	- Turn the device off and on again - check C7
iA	Multifunction input alarm	automatic	check i5 and i6
iSd	High pressure alarm	manual	- Turn the device off and on again - check i5, i6, i8, i9
LP	Low pressure alarm	automatic	check i5 and i6
C1t	Compressor thermal switch alarm	automatic	check i5 and i6
C2t	Second compressor thermal switch alarm	automatic	check i5 and i6
dFd	End defrost time alarm	manual	- Press a key - check d2, d3 and d11

	N.	PAR.	DEF.	SETPOINT	MIN... MAX.
	1	SP	4.0	setpoint	r1... r2
	N.	PAR.	DEF.	ANALOGUE INPUTS	MIN... MAX.
	2	CA1	0,0	Chamber probe calibration	-25... 25 °C/°F if P4 = 3, air probe calibration
	3	CA2	0.0	Evaporator probe calibration	-25... 25 °C/°F
	4	CA3	0.0	Auxiliary probe calibration	-25... 25 °C/°F
	5	P0	1	Probe type	0 = PTC 1 = NTC
	6	P1	0	Decimals	0 = no 1 = yes
	7	P2	0	Temperature measurement units	0 = °C 1 = °F
	8	P3	1	Evaporator probe function	0 = disabled 1 = defrost + fan 2 = fan

	9	P4	0	Configurable input function	0 = digital input 1 = condenser probe 2 = critical temperature probe 3 = air outlet probe if P4 = 3, temperature regulation = product temperature (CPT)
	10	P5	0	Value displayed	0 = Regulation temperature 1 = setpoint 2 = Evaporator temperature 3 = Auxiliary temperature 4 = Air temperature
	11	P7	50	Inlet air weight for calculated product temperature (CPT)	0... 100 % CPT = $\{[(P7 \times (\text{air inlet T})) + [(100 - P7) \times (\text{air outlet T})] : 100\}$
	12	P8	5	Screen refresh time	0... 250 s : 10
	N.	PAR.	DEF.	REGULATION	MIN... MAX.
	13	r0	3.0	Differential setpoint	1... 15 °C/°F
	14	r1	-2.0	Minimum setpoint	99 °C/°F... r2
	15	r2	15.0	Maximum setpoint	r1... 199 °C/°F
	16	r3	0	Setpoint lock activation	0 = no 1 = yes
	17	r4	0.0	Setpoint calibration in energy saving mode	0... 99 °C/°F
	18	r5	0	Regulation for heating or cooling	0 = for cold 1 = for heat
	19	r6	0.0	Setpoint calibration in over refrigeration	0... 99 °C/°F
	20	r7	0	Duration over refrigeration	0... 240 min
	21	r12	0	Differential positioned in r0	0 = asymmetrical 1 = symmetrical
	N.	PAR.	DEF.	COMPRESSOR	MIN... MAX.
	22	C0	1	Compressor delay after start-up	0... 240 min
	23	C1	0	Delay between the start-up of two compressors	0... 240 min
	24	C2	3	Minimum compressor off time	0... 240 min
	25	C3	0	Minimum compressor on time	0... 240 s
	26	C4	5	Compressor off time during probe alarm	0... 240 min
	27	C5	7	Compressor on time during probe alarm	0... 240 min
	28	C6	80.0	High condensation warning limit	0... 199 °C/°F differential = 2 °C/4 °F
	29	C7	90.0	High condensation alarm limit	0... 199 °C/°F
	30	C8	1	High condensation alarm delay	0... 15 min
	31	C10	0	Compressor operating hours	0... 999 h x 100 0 = enabled
32	C11	10	Second compressor activation delay	0... 240 s	
	N.	PAR.	DEF.	DEFROST	MIN... MAX.

33	d0	8	Automatic defrost interval	0... 99 h 0 = manual only if d8 = 3, maximum interval
34	d1	0	Type of defrost	0 = electrical 1 = hot gas 2 = compressor shut down
35	d2	6.0	Defrost final limit	-99... 99 °C/°F
36	d3	30	Defrost duration	0... 99 min if P3 = 1, maximum duration
37	d4	0	Enable defrost upon start-up	0 = no 1 = yes
38	d5	0	Defrost delay after start-up	0... 99 min
39	d6	2	Value displayed during defrost	0 = Regulation temperature 1 = Screen lock 2 = dEF label
40	d7	1	Dripping time	0... 15 min
41	d8	0	Defrost interval counting mode	0 = Hours in device 1 = Hours in compressor 2 = Hours evaporator temperature < d9 3 = adapted 4 = real time
42	d9	0.0	Evaporation limit for the automatic defrost meter	-99... 99 °C/°F
43	d11	0	Activation of final defrost time alarm	0 = no 1 = yes
44	d15	0	Consecutive ON time for hot gas defrost	0... 99 min
45	d16	0	Pre-dripping time for hot gas defrost	0... 99 min
46	d18	40	Adjusted defrost interval	0... 999 min If compressor ON + evaporator temperature < d22 0 = manual only
47	d19	3.0	Limit for Ds adapted (relative to the optimum evaporation temperature)	0... 40 °C/°F Optimum evaporating temperature - d19
48	d20	180	Consecutive time for defrost with compressor ON	0... 999 min 0 = disabled
49	d21	200	Consecutive time for defrost with compressor ON after start-up and over refrigeration	0... 500 min if (temperature regulation - setpoint) > 10 °C/20 °F 0 = disabled
50	d22	-2.0	Evaporation limit for defrost intervals Meter (relative to optimum evaporating temperature)	-10... 10 °C/°F Optimum evaporating temperature + d22
51	d25	0	Air probe enabled for defrost during evaporator probe alarm	0 = no 1 = yes

	52	d26	6	Defrost interval during evaporator probe alarm	0... 99 h 0 = manual only if d25 = 1
	N.	PAR.	DEF.	ALARMS	MIN... MAX.
	53	A0	0	Select value for high/low temperature alarm	0 = temperature regulation 1 = evaporator temperature
	54	A1	-10.0	Low temperature alarm limit	-99... 99 °C/°F
	55	A2	1	Low temperature alarm type	0 = disabled 1 = relating to the setpoint 2 = absolute
	56	A4	10.0	High temperature alarm limit	-99... 99 °C/°F
	57	A5	1	High temperature alarm type	0 = Regulation temperature 1 = Evaporator temperature 2 = Auxiliary temperature
	58	A6	90	High temperature alarm delay after start-up	0... 240 min
	59	A7	1	High/low temperature alarm delay	0... 240 min
	60	A8	15	High temperature alarm delay after defrost	0... 240 min
	61	A9	15	High temperature alarm delay after closing door	0... 240 min
	62	A10	0	Power failure alarm duration recording (not available on EVJ203 and EVJ204)	0... 240 min
	63	A11	1.0	High and low temperature alarm differential restart	1... 15 °C/°F
	64	A12	0	Alarm type power failure notification (not available on EVJ203 and EVJ204)	0 = HACCP LED 1 = HACCP LED + PF label + buzzer 2 = HACCP LED + PF label + buzzer (if duration > A10)
	N.	PAR.	DEF.	FANS	MIN... MAX.
	65	F0	7	Evaporator fan mode during normal operation	0 = OFF 1 = ON 2 = ON if compressor ON 3 = thermoregulated (with temperature regulation + F1) 4 = thermoregulated (with temperature regulation + F1) if compressor ON 5 = according to F6 6 = thermoregulated (with F1) 7 = thermoregulated (with F1) if compressor ON
	66	F1	10.0	Evaporator fan operation limit	-99... 99 °C/°F
	67	F2	1	Evaporator fan mode during	0 = OFF 1 = ON

				defrost and dripping	2 = according to F0
	68	F3	1	Maximum time with fan OFF	0... 15 min
	69	F4	30	Evaporator OFF time during power save mode	0... 240 s x 10 if F0 ≠ 5
	70	F5	30	Evaporator ON time during power save mode	0... 240 s x 10 if F0 ≠ 5
	71	F6	0	High/low humidity operation	0 = low humidity (with F17 and F18 if compressor OFF, ON if compressor ON) 1 = high humidity (ON)
	72	F7	5.0	Evaporator fan limit after dripping (relative to setpoint)	-99... 99 °C/°F setpoint + F7
	73	F8	2.0	Evaporator fan operation differential limit	1... 15 °C/°F
	74	F9	10	Evaporator fans OFF delay after compressor OFF	0... 240 s if F0 = 2 or 5
	75	F10	1	Condenser fan mode	0 = thermoregulated (with F11) 1 = thermoregulated (with F11) if OFF, ON if compressor ON 2 = thermoregulated (with F11) if compressor OFF, ON if compressor ON, OFF during defrost, pre-dripping and dripping
	76	F11	15.0	Limit for condenser fan ON	0... 99 °C/°F differential = 2 °C/4 °F
	77	F12	240	Condenser fan OFF delay after compressor OFF	0... 240 s if P4 ≠ 1
	78	F17	60	Evaporator fan OFF time with low humidity	0... 240 s
	79	F18	10	Evaporator fan ON time with low humidity	0... 240 s
	N.	PAR.	DEF.	DIGITAL INPUTS	MIN... MAX.
	80	i0	5	Door switch input function	0 = enabled 1 = compressor + evaporator fan OFF 2 = evaporator fan OFF 3 = interior light ON 4 = compressor + evaporator fan OFF, interior light ON 5 = evaporator fan OFF + interior light ON
	81	i1	1	Door switch input activation	0 = with closed contact 1 = with open contact

	82	i2	5	Door opened alarm delay	-1... 120 min -1 = disabled
	83	i3	15	Maximum time inhibition regulation with open door	-1... 120 min -1 = until closing
	84	i4	0	Enable alarm recording Door opened	0 = no 1 = yes if i2 ≠ -1 and then i2
	85	i5	8	Multiple uses input function	0 = disabled 1 = power saving 2 = iA alarm 3 = iSd alarm 4 = loaded oper. key 1 ON 5 = loaded oper. key 2 ON 6 = device ON/OFF 7 = LP alarm 8 = C1t alarm 9 = C2t alarm
	86	i6	1	Multiple uses input activation	0 = with closed contact 1 = with open contact
	87	i7	0	Multiple uses input alarm delay	0... 120 min if i5 = 3 or 7, compressor ON alarm delay after restart
	88	i8	0	Number of multiple uses input activations for high temperature alarm applications	0... 15 0 = disabled if i5 = 3
	89	i9	240	High pressure alarm meter reset	1... 999 min
	90	i10	0	Time for consecutive door closings for energy saving mode	0... 999 min After regulation temperature < SP 0 = disabled
	91	i13	180	Number of door openings for defrost	0... 240 0 = disabled
	92	i14	32	Consecutive door opening time for defrost	0... 240 min 0 = disabled
	N.	PAR.	DEF.	DIGITAL OUTPUTS	MIN... MAX.
	93	u1c	0	K1 Output configuration	0 = compressor 1 1 = compressor 2 2 = evaporator fans 3 = condenser fans 4 = defrost 5 = inside light 6 = anti-fog 7 = door resistance 8 = neutral zone resistance 9 = loaded oper.-key 1 ON 10 = loaded oper.-key 2 ON

					11 =alarm 12 =ON/stand-by
94	u2c	12	K2 output configuration		0 = compressor 1 1 = compressor 2 2 = evaporator fans 3 = condenser fans 4 = defrost 5 = inside light 6 = anti-fog 7 = door resistance 8 = neutral zone resistance 9 = loaded oper.-key 1 ON 10 = loaded oper.-key 2 ON 11 =ON/stand-by 12 =alarm
95	u3c	5	K3 Output configuration		0 = compressor 1 1 = compressor 2 2 = evaporator fans 3 = condenser fans 4 = defrost 5 = inside light 6 = anti-fog 7 = door resistance 8 = neutral zone resistance 9 = loaded oper.-key 1 ON 10 = loaded oper.-key 2 ON 11 =alarm 12 =ON/stand-by
96	u4c	2	K4 output configuration		0 = compressor 1 1 = compressor 2 2 = evaporator fans 3 = condenser fans 4 = defrost 5 = inside light 6 = anti-fog 7 = door resistance 8 = neutral zone resistance 9 = loaded oper.-key 1 ON 10 = loaded oper.-key 2 ON 11 =alarm 12 =ON/stand-by
99	u2	0	Enable interior light and operating button on stand-by		0 = no 1 = yes manual
100	u4	1	Enable alarm output buzzer mute		0 = no 1 = yes
101	u5	-1.0	Limit for Door ON Resistance		-99... 99 °C/°F differential = 2 °C/4 °F
102	u6	5	Anti-fog duration ON		1... 120 min

	103	u7	-5.0	Limit for Neutral Zone Resistance (relative to setpoint)	-99... 99 °C/°F differential = 2 °C/4 °F setpoint + u7
	104	u9	1	Enable alarm buzzer	0 = no 1 = yes
	N.	PAR.	DEF.	REAL-TIME CLOCK	MIN... MAX.
	105	Hr0	0	Enable Clock	0 = no 1 = yes
	N.	PAR.	DEF.	ENERGY SAVING	MIN... MAX.
	106	HE2	0	Maximum duration of power saving mode	0... 999 min -1 = until the door closes
	N.	PAR.	DEF.	REAL TIME ENERGY SAVING	MIN... MAX.
	107	H01	0	Energy saving time	0... 23 h
	108	H02	0	Maximum recording time energy saving	0... 24 h
	N.	PAR.	DEF.	REAL DEFROST TIME (if d8 = 4)	MIN... MAX.
	109	Hd1	h-	1st daily defrost time	h- = enabled
	110	Hd2	h-	2nd daily defrost time	h- = enabled
	111	Hd3	h-	3rd daily defrost time	h- = enabled
	112	Hd4	h-	4th daily defrost time	h- = enabled
	113	Hd5	h-	5th daily defrost time	h- = enabled
	114	Hd6	h-	6th daily defrost time	h- = enabled
	N.	PAR.	DEF.	SAFETY	MIN... MAX.
	121	POF	1	Enable ON/STAND-BY key	0 = no 1 = yes
	122	Loc	1	Enable Keyboard Lock	0 = no 1 = yes
	123	PAS	-19	Password	-99... 999
	124	PA1	426	Level 1 password	-99... 999
	125	PA2	824	Level 2 password	-99... 999

9. COOLING MATERIALS

This appliance can be used for chilling all types of materials. The grid shelves can be moved to accommodate different height packages or objects.

To remove the shelves the door must be opened wider than 90°. Substances which give off or absorb odours and flavours as well as liquids should always be stored covered or in closed containers. High-percentage alcohols should be tightly sealed and stored upright.

The ventilation slots on appliances with recirculated air fans on the inside must be not covered.

10. SHUTTING YOUR APPLIANCE DOWN

If your appliance is to be shut down for some time, switch it off and disconnect the plug. Clean the appliance and leave the door open in order to prevent unpleasant smells.

11. DEFROSTING

The refrigerator compartment defrosts automatically. The water that forms on the rearwall drains into a reservoir at the back of the appliance and evaporates automatically through the compressor heat. Ensure that the water can flow without hindrance through the drain hole (see "Cleaning").

A layer of frost or ice can form on the cooling elements when the appliance has been operating at the coldest setting for a long time.

If a **layer of ice** has formed, the cooling elements should be defrosted and the inner surfaces cleaned.

Do not use any mechanical devices or other artificial aids for defrosting other than those recommended by the manufacturer.

12. CLEANING

Periodic cleaning of the device

To clean the device, follow the instructions below:

- Disconnect the equipment from the mains and remove all products stored in it.
- Open all doors and allow time for the interior to reach room temperature. Remove all interior accessories and clean them with soap and lukewarm water. Dry all accessories completely with a soft cloth.
- Once the chamber has reached room temperature, clean all interior and exterior surfaces with soapy water.

Rinse thoroughly and dry with a soft cloth. Failure to dry the device properly may result in water stains.

Stainless steel cleaners are also available in the market, which can repair and protect the protective layer of steel surfaces.

- Put the accessories back in place and connect the unit to the mains.
- Pitting corrosion or cracks in steel are signs of material deterioration. In this case, apply stainless steel cleaners capable of repairing the steel passivation.

WARNING: Never use steel scrubbers, wire brushed or spatulas to clean the device.

NOTE: Cleaning products used must be alkaline-based or chlorine-free. Any cleaner containing chlorides will damage the protective layer of the stainless steel.

Rubber gaskets maintenance

- Rubber gaskets require regular cleaning to protect their elasticity, to ensure proper sealing, and prevent mould growth. Rubber gaskets can be cleaned with soapy water. Avoid using abrasive cleaners or sharp utensils.
- Rubber gaskets can be easily removed in case replacement is needed by pressing them against the door frame.

Cleaning of the condenser

The condenser, located behind the back grille of the device, should be checked regularly. Cleaning frequency will depend on the working environment. Air must flow freely through the condenser, so its surface must be free of dirt and grease. Dirty condensers cause compressor fault and product loss. If the condenser battery is dirty or blocked, follow the steps below:

- Disconnect the device from the mains.
- Remove the back grille from the unit.
- On some models, it will be necessary to remove the screws that affix the condensing unit to the skirting board and remove it in order to clean the condenser.
- If the condenser has a protective case, it must be unscrewed and removed.
- Once the surface of the condenser is accessible, it should be cleaned using a vacuum cleaner or a soft brush. Do not use metal brushes.
- If the dirt is excessive, compressed air can be used for cleaning.
- After cleaning, reinsert the protective case, return the condensing unit to its original position, and re-attach all screws.
- Finally, replace the rear grille and connect the device to the mains.

WARNING!: Do not use water to clean the condenser, as it may damage nearby electrical components.

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Over time and due to normal wear, door hinges may move off place slightly. If you notice that the door is not properly aligned, tightened the screws that attach the hinge brackets to the cabinet.

13. TROUBLESHOOTING

Many operational issues derive from causes that can be easily eliminated without the need to contact the Technical Support Service. The following list covers several types of issues and how to solve them.

PROBLEM	POSSIBLE CAUSE	POSSIBLE SOLUTION
The device does not work	<ol style="list-style-type: none"> 1. The plug is not connected to the power outlet. 2. The plug is not carried with power because the fuse has blown or the automatic current limiter has tripped. 	<ol style="list-style-type: none"> 1. Connect the plug to the power outlet and check that there is electrical power. 2. Replace the fuse or reconnect the automatic power limiter.
The device does not cool sufficiently	<ol style="list-style-type: none"> 1. Check the setpoint temperature on the controller. 2. Door not properly closed or frequent openings. 3. Obstruction of the ventilation grilles of the device. 4. Condenser dirty. 5. The device is exposed to direct sunlight or another heat source. 	<ol style="list-style-type: none"> 1. Decrease the cutting temperature. 2. Take care that the door does not remain open for a long time. 3. Keep these areas free according to the "installation" chapter of this manual. 4. Clean with compressed air or a hard hair brush (not steel). 5. Change the location of the refrigerator or protect it from these heat sources.
Noisy operation	<ol style="list-style-type: none"> 1. The device has not been correctly levelled. 2. Some of the tubes inside the device are rubbing against each other. 3. Loose screws on any of the parts. 4. Fan in condenser or evaporator causing vibrations. 5. Oil load in compressor too low. 6. Loose parts in the condensing unit. 	<ol style="list-style-type: none"> 1. Level as indicated in the "installation" chapter of this manual. 2. Separate the rubbing tubes. 3. Tighten the loose screws. 4. Level the equipment and tighten any loose screws. 5. If the equipment has been overturned at any time, leave it for 24 hours in vertical position without connecting so that the oil returns to the compressor. Check for oil leaks.
The device creates excessive ice in the evaporator	<ol style="list-style-type: none"> 1. Doors not properly closed. 2. Excessive door openings. 3. Defrost has not been carried out. 	<ol style="list-style-type: none"> 1,2. Take care that the door does not remain open for a long time. 3. Schedule a defrost.
Compressor does not start	<ol style="list-style-type: none"> 1. Opened switch. 2. Blown fuse. 3. Faulty wiring. 4. Opened Clixon. 5. Opened controller contacts (defective controller, or device located in a too cold area). 6. Defective relay. 7. Low gas load in the system. 	<ol style="list-style-type: none"> 1. Close the switch. 2. Replace the blown fuse. 3. Check electrical wiring. 4. Check if there is an abnormally low voltage at the socket. 5. Defective controller, or device located in an area that is too cold. 6. Replace the relay. 7. Check for leaks.

Compressor starts, but it shuts down due to overload	<ol style="list-style-type: none"> 1. Low voltage. 2. Faulty wiring. 3. Defective start capacitor. 4. Defective compressor. 	<ol style="list-style-type: none"> 1. Check if there is an abnormally low voltage in the socket 2. Check wiring and electrical installation of the equipment. 3. Replace the start capacitor. 4. Replace the compressor.
High Condensation Pressure	<ol style="list-style-type: none"> 1. Overloaded unit. 2. Air or non-condensable gases in the system. 3. Condenser dirty. 4. Defective condenser fan. 5. Device located in a too hot area. 6. Obstruction in expansion valve or filter. 7. Discharge valve partially closed. 8. Discharge line obstruction. 	<ol style="list-style-type: none"> 1. Check the product temperature, and let it cool outside the equipment if it is very hot. 2. Make a vacuum to extract air or gases. 3. Clean the condenser. 4. Replace the fan. 5. Remove equipment from very hot areas. 6. Adjust expansion valve or filter. 7. Adjust discharge valve. 8. Review and adjust download line.
Reduced Condensation Pressure	<ol style="list-style-type: none"> 1. Insufficient refrigerant load. 2. Leaks in the system. 3. Device located in a too cold area. 	<ol style="list-style-type: none"> 1. Check for refrigerant gas losses. 2. Repair leaks in the system. 3. Remove the equipment from very cold areas.
The compressor performs short cycles	<ol style="list-style-type: none"> 1. Differential control set at too short intervals. 2. Low refrigerant load, check the pressure. 3. Excessive refrigerant load. 4. Leaks in the discharge valve. 5. Opened high pressure switch. 6. Condenser dirty. 	<ol style="list-style-type: none"> 1. Adjust differential with the controller. 2. Check the refrigerant pressure. 3. Adjust the coolant level. 4. Change the discharge valve. 5. Adjust refrigerant charge to avoid excess overpressure. 6. Clean the condenser.