**7000385**

**Agitador de varilla SE-100  
SE-100 rod stirrer**

# Indice

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## 1 Safety

The equipment includes appropriate security measures.

Risky situations that must be respected are indicated along this manual.

### 1.1 Safety icons

The safety icons identify the risky situations and the safety measures that must be taken.

Risk of danger  
Electrical risk

Respect the indicated instructions to perform the operation.

### Electrical risks

Risk of electrical accident when accessing the areas indicated by this sign or when performing the operations indicated by this icon in this manual.

Respect the indicated instructions to perform the operation.

### Important information

- Important information for good results or for an optimum equipment operation.
- Important information to extend life of the equipment or to prevent degradation of any of the components.

### 1.2 Personnel qualification

This equipment can only be used by personnel who have been properly qualified.

This equipment can only be used by personnel who have read and understood these instructions or who has been properly qualified in the operation of this equipment.

## 2 General information

Handle the package with care. Unpack and check that the contents coincides with the "Packing List". If any part is damaged or missing, please advise the distributor immediately.

Do not install or use the equipment without previously reading this instructions manual. They are an inseparable part of the equipment and should be available to all users of the machine. Any questions can be clarified by contacting J.P. SELECTA's technical service.

Any modification, deletion or lack of maintenance of the device could cause damages. The manufacturer is not liable for the damages that may arise.

Do not use the equipment with fluids which can give off vapors or form explosive or flammable mixtures.

### Maintenance

**Manipulation of electronic circuits by unauthorized persons may cause damages difficult to repair. Check with J.P. SELECTA's technical service.**

To clean the accessible parts, use a cloth stained with alcohol cleaner or a similar product. Do not use aggressive solvents.

## Warranty

This product is guaranteed for one year. The warranty does not cover damages caused by an improper use or causes beyond J. P. SELECTA, s.a.u.

Any manipulation of the equipment by an unauthorized personnel cancels the warranty benefits.

### 3 Technical specifications

**Voltage supply 230V 50 Hz as indicated in the machine's features plate.**

Code:	7000385
<hr/>	
External dimensions (H x W x D) mm	175 x 75 x 135
Power	100W
Speed range	60 ... 2000 rpm
Stirring maximum viscosity (mPas)	100.000 cPs
Weight	2.8 Kg
Rod diameter (mm Ø)	Up to 8.5 (through axis)
Maximum torque	450Ncm

### 4 Packing list

#### 4.1 Standard package content

The standard equipment consists of the following components:

- |                       |        |
|-----------------------|--------|
| - Stirrer             | 700385 |
| - Instructions manual | 80320  |

### 5 Introduction

Stirrer 7000385 is a desktop rod stirrer for an accurate and reproducible stirring.

### 6 Installation & Operation

- Place the equipment near a power socket and plug it.
- Place the equipment over a flat level surface, strong enough to support the equipment weight.
- Select the appropriate accessories (see section 7).
- Activate the main switch ON/OFF.
- Select the appropriate speed agitation.  
CAUTION:  
It is advisable to start the agitation with the speed knob at minimum to avoid splashing.

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