<b>TERMOSTATO REGULADOR (S/ADAPT) TV-1452</b>	<b>3001450</b>
<b>TERMOSTATO VISCOSIMETROS TV-1452</b>	<b>3001452</b>
<b>SOPORTE COPA SIN REGULADOR Ø58mm</b>	<b>3001453</b>
<b>SOPORTE COPA SIN REGULADOR DIN Ø58MM</b>	<b>3001454</b>
<b><i>REGULATING THERMOSTAT (W/ADAPTER) TV-1452</i></b>	<b><i>3001450</i></b>
<b><i>THERMOSTAT VISCOSIMETERS TV-1452</i></b>	<b><i>3001452</i></b>
<b><i>CUP SUPPORT WITHOUT REGULATOR Ø58mm</i></b>	<b><i>3001453</i></b>
<b><i>CUP SUPPORT WITHOUT DIN REGULATOR Ø58mm</i></b>	<b><i>3001454</i></b>

**INDICE**

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## GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts, please inquiries calling to your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 2009/104/CE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

## PACKING LIST

The standard equipment consist of the following components:

Description	Code	Quantity
Cup support	(depends on the model)	1
Controller	(depends on the model)	1
Wire	---	1
Level disk	1001027	1
Allen wrench	---	1
Instruction Manual	80180	1

## EQUIPMENT DESCRIPTION

The complete equipment (3001452) is composed by a temperature control unit which includes the keyboard and the temperature reading display and a heating and cooling unit with peltier cells designed to control the temperature of ford cup viscometers.

The regulating equipment consists only of the regulating unit and various adapters, including the heating / cooling unit.

## TECHNICAL FEATURES

Voltage supply 115-230 V 50/60 Hz as indicated on the characteristics plate.

Potencia Power	130 W		
Rango de temperatura Temperature range	10...60 °C		
Estabilidad Stability	±0,1°C		
Error de temperatura Temperature error	±0,5°C		
Resolución Resolution	0,1°C		
Dimensiones controlador Alto/Ancho/Fondo (cm) Controller dimensions Height/Width/Depth (cm)	160	240	215
Dimensiones máx. soporte Alto/Ancho/Fondo (cm) Max. Support dimensions Height/Width/Depth (cm)	235	340	300

## INSTALLATION

Room temperature is decisive when we have to measure newtonian liquids viscosity. So it is highly advisable to install the equipment into a room where the temperature could be controlled to get that it does not differ more than 10°C of test temperature.

Place controller and viscosimeter support on an horizontal surface. Connect both elements by the supplied wire. Plug it to the mains.

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**Be sure that the voltage supply is the same as the one indicated on the characteristics plate.**

**Do not use the equipment if it is not earthed.**

**If you have to change the plug cap be in mind the following:**

**Blue cable: Neutral.**


**Brown cable: Phase.**

**Yellow/green cable: Earth.**

## OPERATION

CONTROL PANEL:

1. Main switch with green indicator lamp.
2. "Set" key
3. "Increase value" key
4. "Reduce value" key
5. Display

**REAR PART:**

6. Fuse carriers
2. Mains wire
3. Connection to viscosimeter support socket

**SETTING EQUIPMENT:**

1. Place the viscosimeter into the support (see picture B). To improve thermic contact it is necessary to use silicone paste for semi-conductors (see picture A). Fix the cup by using an allen wrench number 6 (see picture C).
2. Level out the stand. to check whether the level is correct, use the glass level supplied and adjust it by using the extending legs (see picture D).

Picture A. Silicone paste application

Picture B. Placing viscosimeter

Picture C. Cup fixing


Picture D. Level adjustment

### EQUIPMENT USE:

1. To do viscosity measurements, a stopwatch is needed as well as other supplied elements.
2. Warming cup up: Start up controller switching on main switch (1). Select a temperature pushing the SET key (2) and, at same time, (3) or (4) keys until the display shows required temperature.
3. Filling cup: Block viscosimeter bottom hole with a finger. Fill it with the liquid as well as it spill by the rim.
4. Expelling surplus liquid: Using level glass drive top liquid paying attention that any air bubble remain into the cup. Leave glass as a cover lid. In this moment, bottom hole can be unblocked (liquid will not fall).
5. It can be leaved some minutes to sure that liquid is enough tempered.
6. Place a vessel under bottom viscosimeter hole. Remove glass and at same time start up stopwatch. Stop it the first time liquid trickle cuts.

It is important to remove glass making an horizontal movement as you can see on picture E.

Picture E.

## CALCULATIONS

### VISCOSITY AND ACCURACY CALCULATION

If knowing emptying time, it is possible to calculate kinematic viscosity (see table).

These formulas give viscosity in cSt. To convert units, use the following converting factors.

$$\begin{aligned}
 1 \text{ centistokes} &= 0.01 \text{ stokes} \\
 1 \text{St} &= 1 \times 10^{-4} \text{ m}^2/\text{s} \\
 1 \text{m}^2/\text{s} &= 10000 \text{ St} \\
 1 \text{cSt} &= 1 \times 10^{-6} \text{ m}^2/\text{s} = 1 \text{mm}^2/\text{s} \\
 1 \text{ m}^2/\text{s} &= 1\,000\,000 \text{ cSt}
 \end{aligned}$$

Relationship between viscosity and kinematic viscosity:

$$\text{Centipoises} = \text{centistokes} \times \text{density}$$

The accuracy of measurement, if this has been carried out according to the applicable standard, is 3% (10% on ASTM D 1200).

Taking measurements in the right temperature conditions requires a lot of preparation.

#### Calculation table

Cup type	Time margin (s)	Calculation formula	Measurement range (cSt)	Nominal error
DIN 53211 Ø4	25...150	$v = 4.57t - (452/t)$	90...700	3%
ISO 2431 Ø3	30...100	$v = 0.443t - (200/t)$	5...42	3%
ISO 2431 Ø4	30...100	$v = 1.37t - (200/t)$	35...135	3%
ISO 2431 Ø5	30...100	$v = 3.28t - (220/t)$	100...350	3%
ISO 2431 Ø6	30...100	$v = 6.90t - (570/t)$	190...680	3%
ASTM D-1200 Ø2.53	40...100	$v = 1.44(t - 18)$	27...110	10%
ASTM D-1200 Ø3.40	20...100	$v = 2.31(t - 6.58)$	31...215	10%
ASTM D-1200 Ø4.12	20...100	$v = 3.85(t - 4.49)$	60...370	10%

## SPARE PARTS

15480	Fuse carrier
15511	Male connector
15512	Female connector
20100	Switch
24057	Fan
29293	Control circuit board
29301	Power circuit board
43023	Peltier cell
43066	Pt100 probe
45015	Supply filter
45016	Transformer

## MAINTENANCE


If it is necessary to replace the viscosimeter, you should clean external part of the cup by using a cellulose paper.

Never left that measured liquid dry into viscosimeter. To clean it you must use an appropriate solvent to liquid measured. Do not use metallic tools that can scratch inner cup.

## GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

### Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.